

Kemi Incubator



Product Description

These are used for drying and staining of slides, paraffin embedding, tissue culture work, incubators of anti-body test, excellent for microbiological determination etc.

Power Requirements

Supplied with cord plug switch to operates on 220V, single phase 50Hz, Ac supply.

They are available in following models

Model No.	Innersize
KIA-1/KIS-1	12*12*12 Inches
KIA-2/KIS-2	14*14*14 Inches

Kemi Hot Air Oven



Product Description

These are designed to destroy bacteria, viruses, and fungus and sterilize surgical instruments, glass wares etc by application of dry heat with Temperature ranging from 50°C - 250°C with an accuracy of +/- 1°C.

Power requirements Supplied with cord, plug suitable to operates on 220V single phase, 50Hz Ac supply.

They are available in following models

Model No.	Inner Size
KOA-1/KOS-1	12*12*12 Inches
KOA-2/KOS-2	14*14*14 Inches